

OPERATION INSTRUCTIONS  
UPUTE ZA RUKOVANJE  
GEBRUIKSINSTRUCTIES  
INSTRUÇÕES DE UTILIZAÇÃO  
DRIFTSVEJLEDNING  
BRUKSANVISNING

# Swiss LithoClast® Master

**Rx ONLY** Caution! Federal (USA) law restricts this device to sale by or on the order of a physician

**EMS** 

## COMPONENTS

- 1 Control unit
- 2 Pinch valve
- 3 Holder for stone catcher
- 4 Rotary knob for ultrasound lithotripsy power setting
- 5 Ultrasound standby/ON switch with control lamp
- 6 Connector to Ultrasound handpiece (applied part)
- 7 Pulse mode selection [single pulse/ multiple pulse]
- 8 Toggle mode switch for increase/decrease buttons
- 9 Decrease button (in selected mode)
- 10 Increase button (in selected mode)
- 11 Connector to pneumatic handpiece (applied part)
- 12 Rotary knob for pneumatic lithotripsy power setting
- 13 Pneumatic standby/ON switch with control lamp
- 14 Ultrasound mode display screen
  - a Ultrasound power level [%]
  - b Ultrasound duty cycle setting [%]
  - c Ultrasound use since last stand by [sec/hr]
- 15 Pneumatic mode display screen
  - d Pneumatic power level [%]
  - e Number of impulses since last stand by [sec/hr]
  - f Pulse frequency [Hz]

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
- 16 Power switch ON/OFF
- 17 Power socket
- 18 Fuse box
- 19 Connector for the foot pedal
- 20 Connector for compressed air supply
- 21 Manufacturer identification plate
- 22 Dual foot pedal (combined ultrasound/pneumatic)
- 23 Single foot pedal (for ultrasound only)
- 24 Spare fuses
- 25 Power cord
- 26 Compressed air hose
- 27 Roller clamp
- 28 Autoclavable stone catcher

EMS supplies units with various accessories. The packing list shows exactly what your unit includes.

- 16 Tipka za uključivanje/isključivanje
- 17 Mrežna utičnica
- 18 Kutija s osiguračima
- 19 Priključak za nožnu pedal
- 20 Priključak za dotok komprimiranog zraka
- 21 Identifikacijska pločica proizvođača
- 22 Dvostruka nožna pedala (kombinacija ultrazvučna/pneumatska)
- 23 Jednostruka nožna pedala (samo za ultrazvuk)
- 24 Rezervni osigurači
- 25 Kabel napajanja
- 26 Crijevo za komprimirani zrak
- 27 Rotirajuća stezaljka
- 28 Posuda za skupljanje kamenca s mogućnošću sterilizacije u autoklavu

EMS uz jedinice isporučuje razne dodatke. Popis sadržaja pakiranja prikazuje sve što jedinica obuhvaća.

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## DEAR CUSTOMERS,

Thank you for purchasing this new EMS product. It meets the highest quality and safety standards.

The Swiss LithoClast® Master is a combination ultrasonic and pneumatic urological lithotripter used for the fragmentation of urinary tract calculi in the kidney, ureter, and bladder.

### Contraindications:

Use of the product is contraindicated in patients with the following:

- Active bleeding disorders,
- Solitary functioning kidney,
- Creatinine greater than or equal to 3 µg %,
- During pregnancy,
- Under the age of 18
- Stricture and obstruction problems,
- With an electrical stimulator implanted e.g. pacemaker.

### Potential Complications:

Potential complications associated with fragmentation of urinary tract calculi by pneumatic and/or ultrasound energy are:

- Perforation
- Hemorrhage
- Lesion
- Stone migration
- Pain/colic
- Macroscopic hematuria
- Infection
- Ureteral obstruction

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Instructions concerning the specific components (handpieces, probes) are provided separately in the relative instruction for use. Read the instruction manual delivered with each of the components before using the Swiss LithoClast Master. Do not hesitate to contact your local EMS dealer for further information.

Upute koje se odnose na određene komponente (ručne jedinice, sonda) navedene su odvojeno u uputama za korištenje pojedinačne komponente. Prije korištenja uređaja Swiss LithoClast Master pročitajte priručnik koji se isporučuje sa svakom komponentom. Za više informacija slobodno se obratite svom lokalnom distributeru tvrtke EMS.

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## TECHNICAL MODIFICATIONS

EMS reserves the right to modify the technique, accessories, operating instructions or contents of the set due to technical or scientific improvements.

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## PLEASE READ THIS BEFORE BEGINNING WORK!

### About this manual

Note that the English version of this manual is the master from which translations derive. In case of any discrepancy, the binding version is the English text.

These operating instructions are to ensure the correct installation and use of this product. Always keep this instruction close at hand.

 Please read these operating instructions carefully as they explain all the most important details and procedures. Please pay special attention to the safety precautions.

 This instruction manual describes combined ultrasonic and pneumatic lithotripsy, as well as ultrasonic or pneumatic single lithotripsy application modes.

If your Swiss LithoClast Master has been purchased with the pneumatic lithotripsy mode deactivated, only the ultrasound lithotripsy mode can be used. Please refer in this case only to the instructions for ultrasonic lithotripsy.

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To prevent injury to people and damage to property, please heed the corresponding directives. They are marked as follows:

**Danger**

Risk of patient or user injury

**Caution**

Risk of damage to the device or environmental harm

**Please note**

Useful additional information and hints

**Prohibited**

**Authorized**

We would be pleased to answer your questions or listen to your suggestions. We do, of course, provide support in case of technical problems. Please call directly your approved EMS repair center or your dealer.

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Warnings .....	40
Post-treatment procedure .....	42
Cleaning, disinfecting and sterilizing .....	44
Preparation and cleaning .....	44
Maintenance.....	48
Fuses replacement .....	50
Safety precautions .....	52
Storing of the product .....	56
Product disposal .....	58
Warranty.....	58
Accessories.....	58
EMS customer Service .....	60
Symbols .....	62
Electromagnetic compatibility.....	68
Technical data .....	70
Troubleshooting .....	76

Upozorenja .....	40
Postupak nakon tretmana .....	42
Čišćenje, dezinfekcija i sterilizacija.....	44
Priprema i čišćenje.....	44
Održavanje.....	48
Zamjena osigurača .....	50
Mjere opreza .....	52
Skladištenje proizvoda.....	56
Odlaganje proizvoda na otpad.....	58
Jamstvo .....	58
Dodaci .....	58
Služba za korisnike tvrtke EMS .....	60
Simboli .....	62
Elektromagnetska kompatibilnost.....	68
Tehnički podaci .....	71
Uklanjanje kvarova.....	77

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
## INSTALLATION

 Clean, disinfect and sterilize all the components before first use and after all subsequent uses as described in the "Reprocessing Instructions" manual. Assemble the application parts under sterile conditions.

 Perform function testing according to the instructions given in the respective accessory operating instructions manual.

 Do not cover the unit. Leave a free space of 5 cm around the device for cooling.

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## CONNECTIONS

Connect the power supply cord to the unit. To prevent damage to the unit, check that its rated voltage meets the local line voltage. Connect only to a FI protected mains supply (FI = Residual current protection).

Connect the quick connector of the compressed air tube to the unit, by pushing the connector until it clicks on.

Use only medical-quality, dry compressed air, filtered at a minimum of 5 µm. This supply of compressed air can be either:

- the central compressed air supply system (wall socket),
- a separate compressor,
- a cylinder (ISO 7396).

Connect the foot pedal cord to the unit. Pay attention to the connector orientation and fasten the securing nut.

 Be sure that the connectors are perfectly dry before making the connection. To disconnect, pull only the connector sleeve; never pull on the cord itself.

 When applicable connect the equipotential plug to other medical devices used.

The equipotential plug is used to prevent differences of potential between medical equipments that could generate measurement or treatment errors, if the medical device generates interferences.

In case of observed interferences, the unit shall be connected to other devices used through this plug.

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
**Disposable stone catcher device  
(provided sterile)**

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The disposable Stone Catcher device is provided sterile. Read the instruction manual delivered with it before performing installation.

Connect the short silicone tubing end with the conical connector to the external vacuum source.

 The pinch valve default state is closed and opens when the ultrasound (left) foot pedal is depressed halfway.

 Install a roller clamp on the suction tubing to enable control of the suction flow independently of the lower pressure.

 The sterile disposable stone catcher device cannot be reprocessed.

Jednokratni uređaj za skupljanje kamenca isporučuje se u sterilnom stanju. Prije instalacije pročitajte priloženi priručnik.

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### Stone catcher device tubing

 The use of other tubing than the one supplied with the stone catcher device can either create an incomplete occlusion of the tubing in the pinch valve, or the tubing can collapse during use. Use only the stone catcher device tubing provided by EMS.

The suction tube (standard suction tubing with 3 m length and 5 mm internal diameter) must be sterile and single use. Use of a transparent tubing allows for observation of fragment transport to visualize a possible obstruction. This type of suction tube is available from EMS or can be provided by the hospital.

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## OPERATING THE UNIT

### Suction control

 Recommendations on suction settings for lithotripsy using EMS handpieces and probes are provided separately.

 Suction flow can be adjusted by using the roller clamp on the suction tubing.

 EMS recommends using the EMS Swiss LithoPump, especially designed for this application or another vacuum source with lower pressure control.

 Excessive suction can impair the endoscopic view, collapse the organ, or damage the mucosa. Always check that the suction tubing is inserted into the pinch valve and that the pinch valve is functioning properly.

 Make sure that suction is always operational. Continuous suction during ultrasound lithotripsy is required to evacuate fragments and liquids as well as to ensure the cooling of the probe and the handpiece.

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 An interruption of the suction must be avoided. It is operator's responsibility to continuously monitor the suction flow to check that it is operational and that the level of lower pressure is correctly set.

 Fragments blocked in the lumen of the probe and the handpiece may lead to loss of suction and heating of the probe. If blockage occurs, stop lithotripsy to remove fragments from the probe and from the handpiece lumen before continuing.

If suction stops because of a vacuum pump failure, replace the vacuum pump or use an alternative vacuum source before continuing lithotripsy.

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## Ultrasound power

The left rotary knob is used to adjust the ultrasound power from 10 % to 100 % (in 10 % increments).

## Duty cycle

The duty cycle is the ratio between ultrasound ON and ultrasound OFF over an interval of time of 80 milliseconds. Set the duty cycle using the arrow keys.

The last setting of the duty cycle is kept in memory at power off.

 The handpiece and probe error detection is possible only with a duty cycle set at 100 % (see Troubleshooting chapter).

## Activation

- Press the ultrasound foot pedal halfway to activate the suction (pinch valve opens).
- Press the foot pedal completely to activate ultrasound handpiece, probe, and suction.

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### **Pulse modes**

Select either single or multiple pulse mode with the "Litho mode" button.

### **Pulse frequencies**

When the pulse mode is set to multiple, the frequency can be selected with the arrow keys. Frequency is adjustable on a 1 Hz to 12 Hz range. The last setting of the pulse mode and frequency is kept in memory at power off.

### **Pneumatic power**

The right rotary knob is used to adjust the pneumatic power in percent from 0 % to 100 % (in 10 % increments).

### **Activation**

Press the pneumatic foot pedal to release the pneumatic pulses.

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


### Priprema

### Vorbereiding

 Be sure that the control unit is off while connecting the handpiece.

 When pneumatic lithotripsy is used individually, the suction unit LithoVac lv3 with suction tubes can be added for a simultaneous evacuation of fragments. Refer to the pneumatic handpiece instruction manual for details.

 Suction is controlled with the suction control manual valve on the LithoVac lv3 handpiece.

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 For preparation and settings, refer to the chapters on ultrasound and pneumatic lithotripsy.

Activate the ultrasound mode as well as the pneumatic mode on the control unit.

The "Select" button toggles between the adjustment of the duty cycle and the setting of the pneumatic pulse frequency. The function in use is indicated by a dot either in the duty cycle display or in the frequency display. The arrow keys are used to adjust the values.

 Press both foot pedals to activate pneumatic and ultrasound lithotripsy mode with suction simultaneously (place the foot over both pedals).

 The individual lithotripsy modes can be selectively controlled by pressing only the related foot pedal (left for ultrasound / right for pneumatic).

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## ENDOSCOPY

 It may happen that the probe tip is hidden in folds of the edematous tissue. Endoscopes with a 5° to 15° viewing angle offer a better visualization of the probe tip in the near range of the endoscope.

 When using ureteroscopes, it is recommended to use a handpiece guidance to obtain a proper probe to endoscope alignment and to limit the extension of the probes beyond the endoscope tip. Guiding adapters can be purchased from EMS for the most common models of ureteroscopes.

 Ensure that the tip of the pneumatic or ultrasound probe does not extend beyond the tip of the endoscope by more than 10 to 20 mm when it is fully inserted in the endoscope working channel.

 Risk of mucosal perforation: When using a 0.8 mm or 1.0 mm pneumatic probe with an ureteroscope, it is recommended to adjust (if needed) the length of the probe to the endoscope length. Contact the EMS after sales service to adjust the probe length to fit a particular endoscope.

 Adjusting the length of the probe and using a guiding adapter will help limit the risk of perforation and reducing energy loss, particularly with thin probes.

 To avoid damage to the endoscope, the probes should not be activated with their tip still inside the working channel of the endoscope.

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## TREATMENT RECOMMENDATIONS

### General

Intracorporeal lithotripsy by use of ultrasound and pneumatic energies, used simultaneously or individually, requires experience and skill in endoscopic technique. A good balance of irrigation through the endoscope versus suction is important for a proper endoscopic view. Excessive and badly timed flushing should be avoided in order not to flush away the stone.

The stone can only be fragmented if the ultrasound probe and/or the pneumatic probe is/are in direct contact with it. Before releasing a pulse or activating the ultrasound, touch the stone with the probe tip and move it slightly to verify contact.

Due to the push back effect, mobile ureteric stones should not be treated from below with the probe positioned in the center of the stone. Best practice is to position the probe tangentially at the stone so that the lithotripsy pulses will not push it back.

The stone should be hit at the corner or at the edges depending on its structure.

It is useful to treat the stone by pushing it slightly against the ureter wall and working with a low pulse rate. Use of a stone migration prevention device (e.g. Dormia basket) may be indicated.

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It is possible to treat impacted stones directly "in situ" using only pneumatic lithotripsy, because pneumatic lithotripsy does not generate heat, therefore no risk of thermal damage to surrounding tissues.

### Recommended settings

 For each lithotripsy method, whether used individually or in combination, it is recommended to start fragmentation at medium power settings of 50 % to 60 %.

The frequency of the pulses should not be higher than 7 Hz in individual mode with a duty cycle set at 100 % at the beginning whether used alone or in combination with pneumatic lithotripsy, and can be increased or decreased gradually.

For combined application, it is recommended to set the pneumatic pulses frequency at 2~3 Hz to help controlling the fragmentation.

 Working with a duty cycle of less than 100 % can be helpful to reduce warm-up of the ultrasound instrument if e.g. the endoscopic procedure requires to reduce suction through the ultrasound instrument.

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### **Selection between individual and combined modes**

The Swiss LithoClast Master offers two modes of lithotripsy: Ultrasound lithotripsy and pneumatic direct impact lithotripsy, known as the Swiss LithoClast principle.

These two lithotripsy modes can operate in combination or individually. The selection of the treatment mode depends on the stone location and on the stone properties.

Stone fragmentation testing has been conducted to verify that simultaneous use of ultrasound and shock wave pneumatic lithotripsy can reduce the time needed for stone fragmentation and evacuation. In situations where it is time consuming and hard to fragment the stones by using ultrasound energy alone, the simultaneous application of the pneumatic probe guided through the ultrasound probe can considerably improve the process. This is particularly true when treating hard stones, such as mono hydrate stones.

The combined use is not indicated when the fragmentation process works fine with the ultrasound probe alone. The flow characteristic of the ultrasound probe is limited when the pneumatic probe is inside. Fragment evacuation is more efficient when the suction channel is kept free.

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## WARNINGS

 Danger of urothelial and mucosal perforation! Due to the tendency of the ultrasound probe to perforate the stone rather than fragmenting it, there is a risk of injuring the mucosa of the kidney and bladder behind the stone. During the entire treatment, the user must keep the probe tips under endoscopic view (the probes tip is extended 10-20 mm beyond the endoscope tip).

 Do not apply pressure with the distal part of the probes against tissues as damage to tissue or perforation could occur. Use particular caution with thin probes.

 Do not touch the ultrasound probe when ultrasound is activated. The ultrasound vibrations of the probes may cause heat and instant burns.

 The user's hand may have a tendency to push forward under the influence of the repeated pneumatic pulses. This tendency must be corrected by the user who will continuously control the relative position of the probe to reduce risks of perforations.

 When using the unit in multiple pulse mode, be sure to check the progress of fragmentation and the probe position on a frequent basis (every 1 to 2 seconds).

Sterile backup ultrasound and pneumatic lithotripsy probes must be available in the event of a probe failure during treatment.

If a probe breaks distally, use sterile grasping forceps to remove probe pieces from the urinary tract.

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## **POST-TREATMENT PROCEDURE**

### **Emptying of suction tubes and stone catcher**

- Remove the silicone tubing of the stone catcher from the pinch valve,
- Suck air in to empty the suction tubes,
- Tilt the stone catcher in order to empty the receptacle by suction.

### **Elimination or conservation of stone catcher content**

- If fragments are not to be kept for analysis, dispose of the stone fragments per hospital protocol.
- If the fragments are to be kept, close the receptacle with the yellow transport closing cap.

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## CLEANING, DISINFECTING AND STERILIZING

 Please conform to the recommendations of the "Reprocessing Instructions" manual delivered with your product regarding procedure of cleaning, disinfecting, packaging and sterilizing of the components.

 Please follow present-day rules about sterilization applicable in your country.

### Preparation and cleaning

#### Unit

 The housing of the unit is not watertight.

Please proceed as follows:

- Turn off the device,
- Disconnect the power supply connector,
- Disconnect the air supply connector,
- Disconnect the foot pedal connector,
- Disconnect the handpiece(s) and attach the caps onto the connectors.

 Clean the unit only with an alcohol-based, commercially available (ethanol, isopropanol), colorless disinfectant. The use of scouring powder or an abrasive sponge will damage its surface.

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## MAINTENANCE

 Should legal provisions in your country demand maintenance intervals, these must be observed.

### Display of the software version

When the device is turned on, the installed software version is displayed for few seconds.

### Display of cumulated operation time

During power up, the total time of ultrasound usage as well as the total number of pneumatic impulses are displayed in the related displays on front panel.

### Ultrasound time

The cumulative ultrasound usage time [seconds/ hours] is shown briefly on the left display at start-up.

Remark: If the number is superior to 9999 (seconds), the hours mode is shown, for example: 012h = 12 hours. If the number is superior to 999h (hours), the exponential mode is used:  
i.e.  $45E3 = 45 \times 10^3 = 45'000$  hours

### Pneumatic pulse counter

The cumulative number of pulses is shown briefly on the right display at start-up.

Remark: If the number is superior to 9999, the exponential mode is used  
For example:  $56E4 = 56 \times 10^4 = 560'000$

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Test kontrollenheten. Hvis sikringen går igjen, kontakter du forhandleren.

## SAFETY PRECAUTIONS

EMS and the distributor of this product accept no liability for direct or consequential injury or damage resulting from improper use, arising in particular through non-observance of the operating instructions, or improper preparation and maintenance.

 Read this instruction manual before using this device. This product may only be applied for its intended use by qualified personnel and for the applications described in this manual. If the unit is used in combination with other instruments, please refer to their instruction.

 Use for the intended purpose only. Before using the product, make sure the operating instructions have been carefully studied. This also applies to any equipment used with this product. Failure to observe the operating instructions may result in the patient or user suffering serious injury or the product being damaged.

 Risk of explosion: Do not use this product in the presence of inflammable anaesthetics or gases.

 No specific training other than initial professional training is required to use this medical device. The practitioner is responsible for performing the clinical treatments and for any dangers that may arise due to a lack of skill and/or training.

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 Always examine the product for damage before commencing treatment. Damaged accessories or a damaged unit must not be used and must be replaced. Use original EMS spare parts and accessories only.

 Proper care and maintenance are required to ensure safe operation of medical devices. EMS therefore recommends that devices and instruments be checked regularly for proper functioning and completeness, upon receipt and prior to each application, to prevent injury to patients and users.

 Product repairs are to be executed exclusively through an approved EMS repair center. After any repair or modification, the technical safety of the device and accessories must meet the corresponding standard applicable for such equipment. Each repair, modification or test must be recorded in a report.

 Before each use, always clean, disinfect and sterilize the different pieces and accessories of this product. Please refer to the information provided in the "Reprocessing Instructions" manual delivered with your product. Non-sterile pieces and accessories may cause bacterial or viral infections.

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
## STORING OF THE PRODUCT

## SKLADIŠTENJE PROIZVODA

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 Keep the original packaging until the product is to be disposed of permanently. You can use it for shipping or storing at any time.

Should you wish to store the product for an extended period of time:

- Proceed as described in the chapter "Cleaning, disinfecting and sterilizing",
- Pack the product and all accessories in the original packaging.

Please refer to the storage and transport conditions in the "Technical data" section.

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## PRODUCT DISPOSAL

The product must not be discarded in domestic household waste.

Should you wish to definitively dispose the product, please comply with the regulations which apply in your country.

Waste Electrical and Electronic Equipment belonging to customers located in the European Union may be shipped to EMS for recycling in accordance to the WEEE regulations. The costs of recycling, exclusive of shipping fees, are covered by EMS.

## WARRANTY

Damages due to non-adherence to the operating instructions or wear out of parts are excluded from warranty.

 The warranty of your product will be cancelled if you try to open it.

## ACCESSORIES

Accessories are available from EMS or any authorized dealers. Please contact your customer service directly.

 For safety reasons this product may only be used with genuine EMS accessories. The user assumes all risks of using this product with non-genuine accessories or otherwise in non-compliance with the instructions manual.

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## EMS CUSTOMER SERVICE

Please contact your dealer or EMS repair center for any servicing or repairs.

EMS decline responsibility for the safety of the device and declare the warranty null and void if service or repair is carried out by an unauthorized third party or if non-genuine spare parts are used.

It is best to ship your product in the original packaging. It protects your product against damage during shipment.

 Before dispatching your product including all accessories, please clean, disinfect and sterilize as described in the operating instructions.

 Risk of transmitting micro-organisms! To protect the personnel of the repair center and for safety reasons during transport and shipment, all products and accessories returned to the factory for repair or revision must be cleaned, disinfected and sterilized in accordance with the instructions manual. Repair can be refused for product or accessories received in a contaminated condition.

When sending your product directly to the approved EMS repair center, please include the name of the dealer to simplify processing.

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Molimo obratite se svom dobavljaču ili centru za popravke tvrtke EMS za bilo kakvo servisiranje ili popravke.

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## ELECTROMAGNETIC COMPATIBILITY

 This product needs special precautions regarding electromagnetic compatibility and needs to be installed and put into service according to the note "Electromagnetic compatibility".

 Certain types of mobile telecommunication equipment could potentially interfere with this product. The separation distances recommended in the note "Electromagnetic compatibility" must be taken into account.

 This product should not be used adjacent to or stacked with an other unit. If adjacent or stacked use is necessary, functioning should be verified prior to operation in the configuration in which it will be used.

 The use of accessories and cords other than those specified or sold by EMS as replacement parts may result in increased emissions or decreased immunity of this product.

 The equipotential plug is used to prevent differences of potential between medical equipments that could generate measurement or treatment errors, if the medical device generates interferences. In case of observed interferences, the unit shall be connected to other devices used through this plug.

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## TECHNICAL DATA

### DESCRIPTION

Manufacturer	EMS SA, CH – 1260 Nyon, Switzerland
Model	Swiss LithoClast® Master
Power supply (voltage)	100 – 240 VAC
Power supply (frequency)	50 – 60 Hz
Power consumption	500 VA
Ultrasound frequency	23 – 26.4 kHz
Ultrasound power consumption	78 W
Energy output	Approx. 85 mJ
Classification	EN-60601-1: Class I Applied part BF IP 20 (Foot pedal: IP X8) MDD 93/42 EEC: Class IIb
Operating mode	Continuous modes
Primary fuse	T4A 250 V (Ø5 X 20mm, time lag, low breaking capacity)
Compressed air supply	3.5 to 6.5 bar (3500 to 6500 hPa)
Housing	Weight: 11 kg Dimensions: height: 135 mm Width: 360 mm Depth: 420 mm
Operating conditions	Temperature +10°C to +30°C Relative humidity 30 % to 75 % Atmospheric pressure 700 hPa to 1060 hPa Altitude max. 3000 m

Storage and transport conditions	Temperature –10°C to +40°C Relative humidity 10 % to 95 % Atmospheric pressure 500 hPa to 1060 hPa
Shelf life / Validated lifetime	Sterile accessories shelf-life: refer to packaging Handpiece: 100 reprocessing cycles Wench lifetime: 3 years, or 300 reprocessing cycles Console: 7 years

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## TROUBLESHOOTING

### Error codes description

The Swiss LithoClast Master provides an error code system in case of malfunction or incorrect handling. The code is made of two blocks of characters like *Er 12*. On the display, the code is on the left for ultrasound, on the right for pneumatic, and is split left-right for general. Codes of the A type disappear with the cause. The B type code must be confirmed following solution pressing and releasing the foot pedal. The ultrasound duty cycle must be at 100% for the error message to be detected.

No.	type	linked to	cause/description	solution
Er 10	A	push button or pedal	Defective foot pedal or foot pedal accidentally pressed while switching the unit on.	Pay attention not to press the foot pedal when switching on the unit. Check that nothing pushes on or impair the normal functioning of the pedal. If failure persists call EMS after sales service.
Er 11			Defective push button or push button accidentally pressed while switching the unit on.	Pay attention not to press a push button when switching on the unit. If failure persists call EMS after sales service.
Er 20	B	ultrasound handpiece	The internal ultrasound power value cannot be obtained. Defective ultrasound probe or handpiece.	Please proceed as follows: 1 Verify that the ultrasound probe is tightened, 2 Press the foot pedal to try again, 3 If the error appears again, replace the probe, 4 Press the foot pedal to try again, 5 If the error appears again, reduce the ultrasound power to 80 %, 6 Press the foot pedal to try again, 7 If the error persists, Replace the ultrasound handpiece, 8 Press the foot pedal to try again, 9 If the error is still present, call EMS after sales service.
Er 21			Defective ultrasound handpiece or cord (i.e. a damaged insulation).	
Er 23			Loss of ultrasound effectiveness, the control unit cannot adjust the resonance frequency. The handpiece is loose or not connected, the ultrasound probe is damaged or untightened.	
Er 24			Ultrasound efficiency is not stable. The resonance frequency is found and lost (similar to code Er 23 but this code does not come up because the ultrasound functions intermittently).	
Er 30	B	ultrasound mode	Internal power supply failure.	Turn the device off, then on again. If failure persists call EMS after sales service.
Er 40	A	pneumatic mode	Internal working pressure is not in the specified range. Input pressure can be too high or too low. Activation of the pneumatic pulses is stopped.	1 Verify and adjust the input pressure (supply side) according to the value specified on the identification plate on the unit back panel. 2 Reduce working pressure using the pneumatic rotary knob. Press the foot pedal to try again. 3 Do not increase pressure beyond 100 % again. If the error is still present, call EMS after sales service.

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UK

**CAMBRIDGESHIRE, UNITED KINGDOM**  
**EMS UK LTD SUITE 2**  
Xenus House, Sandpiper Court  
Phoenix Park - Eaton Socon  
UK - St Neots, Cambridgeshire  
PE19 8EP  
T +44 01480 587260  
E-mail: info@ems-unitedkingdom.com

AU

**SYDNEY, AUSTRALIA**  
**EMS OCEANIA PTY LTD**  
Suite 3, Shed 73  
4E Huntley Street  
AU - Alexandria, NSW 2015  
T. +61 (02) 9313 4392  
E-mail: info@ems-oceania.com

KR

**SEOUL, SOUTH KOREA**  
**EMS ELECTRO MEDICAL SYSTEMS KOREA LTD.**  
13, Beobwon-ro 8-gil, Songpa-gu  
T. +82 70 7807 2875  
E-mail: emskorea@ems-ch.com

PL

**KATOWICE, POLAND**  
**EMS POLAND SP. Z O.O.**  
Al. Rozdzińskiego 188H  
PL - 40-203 Katowice  
T +48 32 4937060  
E-mail: info@ems-poland.com

© Copyright EMS  
FB-411/4  
Rev.D  
ed. 2022/02

**EMS - SWISSQUALITY.COM**

## Signature Manifest

Document Number: FB-411/4

Revision: D

Title: IFU\_Swiss\_LithoClast\_Master

All dates and times are in UTC+01:00.

### CR-11109 - MED - Addition of Life Cycle information to FB-411/4, FB-560/32, FB-569/4, FB-592/4

#### Author and Reviewers

Name/Signature	Title	Date	Meaning/Reason
Maxime Jost (MJOST)		21 Mar 2022, 12:45:07 PM	Approved

#### Final Reviewer

Name/Signature	Title	Date	Meaning/Reason
Suzanne Fassio Hardy (SFA)		21 Mar 2022, 01:45:53 PM	Approved

#### Approval

Name/Signature	Title	Date	Meaning/Reason
Timothée Deblock (TDE)		22 Mar 2022, 08:58:12 AM	Approved

#### Update ERP

Name/Signature	Title	Date	Meaning/Reason
Marie-Christine Lamy-Biehler (MCL)		22 Mar 2022, 02:05:56 PM	Approved

#### Notification

Name/Signature	Title	Date	Meaning/Reason
Generic user Notificator (NOTIFICATOR)		22 Mar 2022, 02:05:57 PM	Email Sent
Maxime Jost (MJOST)		22 Mar 2022, 02:05:57 PM	Email Sent